

NovoCyte

High performance flow cytometer for everyone

Roman Hudec, PhD

Field application specialist

GSM: +421 907 738 115

E-mail: hudec@accela.eu



Doc. Tomáš Mackulák
FCHFT STUBA Bratislava


31. Jan – 1. Feb 2018

Material and methods

- **ACEA NovoCyte 3000 (488, 642, 405 nm lasers), 13 colors**
physical parameters: FSC and SSC
fluorescence: excitation 405nm laser
emission channel VL2 530/30nm AmCyan
VL3 572/28nm Pacific Orange
- **10g of mechanically grinded polystyrol – phenyl-thiophen copolymer with Mw 250 000 g/mol (Mn 90 000g/mol) in ddH₂O was filtrated through 45µm filter.**
- **10ml of heterogenous mixture was used for degradation reaction using 0, 1.1, 7 and 16.8 mg of Fe⁶⁺**
- **Mixture was analyzed by NovoCyte after 20min of the reaction**

Degradation of fluorescent polystyrol

➤ Physical parameters




➤ Increase in particle counts

Degradation of fluorescent polystyrol

➤ Fluorescence

VL2 | 405 | 530/30 | AmCyan




➤ Decrease in fluorescence pos. and increase in neg. particles

Degradation of fluorescent polystyrol

➤ Fluorescence

VL3	405	572/28	Pacific Orange
-----	-----	--------	----------------



➤ Decrease in fluorescence pos. and increase in neg. particles