

# Application Note



## Cell-by-cell Analysis of Ultra High-density HepG2 Culture

### Introduction


For years, many types of high-density culture methods, such as spheroids, fiber scaffolds, and extracellular matrixes, have been proposed for *in vitro* cell-based assays. These culture systems have been recognized to more accurately simulate a cells natural environment than standard monolayer cultures on a flat substrate. Therefore, cells in high-density culture conditions are expected to exhibit responses against chemical treatments that closely resemble responses of tissues *in vivo*.

General homogeneous assay protocols originally developed for monolayer culture can be applicable to these high-density culture systems with minor modifications. However, for microscopy and other image-based assays, there are significant obstacles to overcome when applying conventional image analysis protocols to high-density cultures. The optical architectures of most microscope-coupled research instruments do not capture light based information throughout the entire thickness cell aggregates due to depth-of-field limitations. In addition, image-analysis software optimized for monolayer cell culture are not able to perform cell-by-cell object recognition and resulting quantification.

Here we show a set of examples with an ultra-high density HepG2 (hepatocellular carcinoma) cell culture to explain how the CQ1 can capture clear images from entire thickness of a cell layer three-dimensionally. Additionally, we will show how the CQ1 can analyze cell responses against chemical treatment, on a cell-by-cell basis.

The protocol on this note has a potential applications for analyses of various high-density cultures, including three dimensional cell culture conditions.

### Negative control



### Staurosporine treated




Figure 1. Images of highly dense HepG2 cultures and nuclear contouring. Non-treated (a) or staurosporine ( $10^{-7}$  M for 48 h) treated (c) cells. Three dimensional cell-by-cell nuclear recognition was carried out by using the Spheroid Analysis Algorithm of the CQ1 software (b and d). Objective lens: 20X.


## Experimental procedure

- A 96-well glass bottom microplate was coated with an extra cellular matrix, Matrigel (5 fold dilution with culture medium).
- HepG2 (hepatocellular carcinoma) cells were seeded on the pre-coated plate at a density of  $5 \times 10^4$  cells/well. The plate was incubated for 48 hours to create an over-confluent state.
- Staurosporine was added to experimental wells and the plate was incubated for an additional 48 hours.
- Cells were fixed with a formaldehyde solution then tagged with anti-active caspase-3 and anti-H3Ser10P. Bound primary antibodies were visualized with fluorescently labelled secondary antibodies. Cell nuclei were stained with Draq7 in the presence of RNaseA.
- Cell images were captured and the images analyzed by the CQ1. Graphing and statistical processing was carried out using FCS Express™ 5 Image Cytometry (De Novo Software, Glendale, CA) (optional).

## Results and discussions


An ultra high-density HepG2 cell culture was created and cell-by-cell analysis was carried out by the CQ1 to evaluate cell responses against staurosporine, a hepatotoxic chemical.

- To facilitate dense cell layer formation, the plate was coated with an extracellular matrix. In preliminary tests, compared with culture on a non-coated normal plastic-bottom plate, HepG2 cells plated on the matrix-coated substrate exhibited approximately ten-fold higher sensitivity toward staurosporine toxicity (data not shown).
- Images of fluorescently labelled cells were captured in three-dimensions from the thick cell layer formation (Fig 2). Two molecular markers, H3Ser10P (cell growth) and active caspase-3 (apoptosis) were selected for immunofluorescent labelling.
- Segmentation of individual cell nuclei enabled cell-by-cell characterization in response to staurosporine (Fig 3).
- The CQ1 is a versatile system that allows simultaneous analysis of multiple markers and parameters at the single cell level in high density cell culture based assays, including thick cultures of hepatocytes.



**Figure 2. 3D reconstruction of HepG2 cell images.**

Twenty-one slices (along the Z-axis) of a multicolor image encompassing 50  $\mu\text{m}$  of thickness were reconstructed to create a 3D image. Non-treated (a) or staurosporine treated (b) cells were fluorescently immunostained with anti-H3Ser10P (magenta) and anti-active caspase-3 (green). Cell nuclei were stained with Draq7 (gray). Objective lens: 20X.



**Figure 3. Multi-parametric analysis of two cell markers**

Image analysis data from the CQ1 of non-treated (a) or staurosporine treated (b) cells was exported and further analyzed by scatter plots in FCS Express™ 5 Image Cytometry. The proportions of the growing or dead cell populations were evaluated quantitatively. On each scatter plot,  $1 \times 10^4$  events were plotted.

**YOKOGAWA** 

YOKOGAWA ELECTRIC CORPORATION

Measurement Business HQ, Life Science Center

Kanazawa 2-3 Hokuyodai, Kanazawa-shi, Ishikawa, 920-0177 Japan  
 Tokyo 2-9-32 Nakacho, Musashino-shi, Tokyo, 180-8750 Japan  
 Osaka Breeze Tower, 4-9 Umeda 2, Kita-ku, Osaka, 530-0001 Japan  
 E-mail: CSU\_livecell\_imaging@cs.jp.yokogawa.com

Phone: (81)-76-258-7028 FAX: (81)-76-258-7029  
 Phone: (81)-422-52-5550 FAX: (81)-422-52-7300  
 Phone: (81)-6-6341-1408 FAX: (81)-06-6341-1426  
 URL: <http://www.yokogawa.com/scanner>